

QUICK START GUIDE
950 WRIST SPLINT
SKU: SUP2038 / SUP2039



APPLICATION INSTRUCTIONS

1. Loosen straps and slide your hand into the brace.



2. Fasten the straps from bottom to top.





Distributed by



8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com

011225

Note: All straps can be cut to size.



CARE INSTRUCTIONS

Remove the main splint. Hand wash with a mild detergent in cold water. Hang dry. Do not machine wash or dry.

WARNINGS ⚠

- Consult your doctor before use.
- Stop use if pain occurs or worsens.
- Assistance may be required during the application of this support.
- Do not over tighten, doing so may lead to restricted blood circulation.
- If you experience any pain, swelling, sensation changes, or any unusual reaction while using this product please consult your doctor immediately.
- Do not machine wash or dry the brace.