



# 807 THUMB SPLINT

SUP2040BLK / SUP2041BLK



To see all the FAQ's,  
visit [coretechortho.link/1bm](https://coretechortho.link/1bm)

Distributed by

**vive**  
health

8955 Fontana Del Sol Way  
Naples, FL 34109  
1-800-487-3808  
[www.vivehealth.com](https://www.vivehealth.com)

**coretech**<sup>®</sup>

## QUICK START GUIDE

## APPLICATION INSTRUCTIONS

1. Loosen all straps and carefully insert your hand so that the thumb is positioned into the splint.
2. Tighten the straps from bottom to top.



**NOTE:** All straps can be cut to size.



## CARE

Remove all splints, hand wash with mild detergent.  
Air Dry

## WARNINGS

- Consult your doctor before use.
- Stop if pain worsens.
- If you experience any pain, swelling, sensation changes, or any unusual reaction while using this product please consult your doctor immediately.
- Do not over tighten, doing so may lead restricted blood circulation.