

www.coretechortho.com

coretech[®]
by vive health

OWNER'S MANUAL
397 SUSPENSION SLEEVE
SKU: SUP2067



To see all FAQs in one place, visit, [coretechortho.link/sx4](https://www.coretechortho.com/link/sx4)

397 SUSPENSION SLEEVE

ADJUSTMENT

1. The Coretech Suspension Sleeve is universally sized for the right or left knee. However, it can be trimmed to a specific length for a more comfortable fit if needed (optional).
2. Ensure the sleeve fits properly before altering with scissors.



APPLICATIONS

To Prevent Skin Irritation:

1. Sit comfortably, and wrap the Coretech Suspension Sleeve around the desired area before the knee brace is applied. Suspension sleeves are most commonly applied around the thigh, upper calf, or directly over the knee.



2. Use the free-sized thumb loop to pull and tighten securely for optimum tension.
3. Ensure the sleeve stays comfortable and secure by attaching the loose end to the hook and loop provided.
4. Once your suspension sleeve is secured, apply the knee brace as usual.



To Prevent Brace Migration:

1. While standing, wrap the suspension sleeve around your fitted knee brace.
2. Secure the sleeve to the brace using the removable hook tabs.
3. Tighten with the thumb loop until you achieve appropriate tension.
4. Fasten the sleeve by attaching the top and bottom material.



CARE

- Hand wash in cold water with mild detergent.
- Leave to air dry.
- Do NOT machine wash or dry.
- Do NOT clean, iron, or bleach any part of the item as it can cause damage.

WARNINGS

- Consult your doctor before use.
- Stop use if pain worsens.
- Assistance may be required during the application of this support.
- Do not over tighten, doing so may lead to restricted blood circulation.
- If you experience any pain, swelling, sensation changes, or any unusual reaction while using this product please consult your doctor immediately.
- Further questions regarding the adjustment or duration of suspension sleeve use should be forwarded to your healthcare provider.



Distributed by

vive
health

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com