

www.coretechortho.com

coretech
by vive health

OWNER'S MANUAL HINGED KNEE BRACE

SKU: SUP2072BLK



To see all FAQs in one place, visit, [coretechortho.link/5pt](https://www.coretechortho.com/link/5pt)

HINGED KNEE BRACE

APPLICATION

1. While seated, position the brace around your knee. The round opening should be aligned with your patella.
2. Wrap the bottom strap below your knee.
3. Wrap the top strap above your knee.
4. Secure the brace with the remaining two straps.



5. Walk around to test the fit and adjust as needed.



CARE

We recommend removing the two hinged splints before washing.



Hand wash in cold water with mild detergent. Hang dry.

WARNINGS !

- Consult your doctor before use.
- Do not continue use if pain worsens.
- Do not overtighten.
- Do not machine wash or dry.

coretech

by **vive health**

Distributed by

vive
health

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com