

www.coretechortho.com

coretech[®]
by **vive health**

OWNER'S MANUAL
386 WALKER BOOT SHORT
SKU: SUP3092BLKS



To see all FAQs in one place, visit [coretechortho.link/59d](https://www.coretechortho.com/link/59d)

386 WALKER BOOT SHORT

PUTTING THE WALKER ON



For a video demonstration,
check out coretechortho.link/074

1. Unfasten all secure straps and open the padded liner.



2. (Optional) Two extra pads are included to cushion the inside of the liner,



they can be used around the ankle, the calf, under the foot, or the straps.

3. With your foot inside the liner, secure the bottom, then the top strap.

Make sure your heel is positioned completely back.



SUP3092BLKS

4. Secure the outer straps from bottom to top and tighten to a comfortable fit.



Do NOT overtighten.



CARE INSTRUCTIONS

- Remove the padded liner. Hand wash in cold water with a mild detergent. Air dry.
- Wipe the walker with a damp cloth and dry immediately.
- Do NOT machine wash or dry.

NOTES



Distributed by

vive
health

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com