

HAVE MORE QUESTIONS?

Check out our list of Frequently Asked Questions at [vhealth.link/n6j](https://vivehealth.link/n6j) for helpful answers.



And if that doesn't answer your question, our customer service team would love to help! Feel free to connect with them by phone, e-mail, or chat on our website.



service@vivehealth.com



1-800-487-3808



vivehealth.com



Carbon Fiber Cane
MOB1039BLK



(01)00810041983343(241)MOB1039BLK

Distributed by

vive
health

8955 Fontana Del Sol Way
Naples, FL 34109
Made in China



[vivehealthusa](https://vivehealthusa.com)

vive®

CARBON FIBER CANE

Owner's Manual
MOB1039

vivehealth.com



To see all FAQs
visit [vhealth.link/n6j](https://vivehealth.link/n6j)

FIND YOUR CORRECT CANE HEIGHT

1. Let your arm hang naturally by your waist. Adjust the cane height so the handle touches your wrist.



WARNING

- Make sure the brass pin is locked into place before use.
- Weight Capacity: 220 LBS

CANE HEIGHT ADJUSTMENT

1. Loosen the locking latch.
2. Slide extension tube up or down to desired height. Button snap should secure cane into place.
3. Tighten the locking latch.

