

# EXTENDED TRIGGER FINGER

SUP1071

## WHAT'S INCLUDED

- 1x Full Finger Splint
- 1x Long Palm Strap
- 1x Short Finger Strap
- 1x Wrist Strap



Check out this how-to video at  
[vhealth.link/8ec14](https://vhealth.link/8ec14)

**vive**  
vivehealth.com

## QUICK START GUIDE

## USING THE SPLINT



For video demonstration  
visit [vhealth.link/501b2](https://vhealth.link/501b2)

1. Place the hand over the splint so that the injured finger lays over the splint.
2. Secure the top strap at the fingertip.



3. Slip on the wrist strap through the top.



Secure via the hook and loop tab on the inner side, tighten the strap through the d-ring.



4. Cut the excess material from the splint end.



QUICK START GUIDE

vive  
vivehealth.com

Distributed by

vive  
health

8955 Fontana Del Sol Way  
Naples, FL 34109  
1-800-487-3808  
www.vivehealth.com

170424

5. Attach and secure the palm strap.



6. Attach and secure the mid-finger strap.



7. Adjust the straps to comfort.



## CARE INSTRUCTIONS

- Hand wash in cold water with a mild soap.
- Rinse thoroughly
- Air dry
- Keep the splint dry while in use.
- DO NOT machine wash or dry.
- The fabric is latex-free and the splint inside is aluminum to resist corrosion.